

EMFIT MM
Movement Monitor™

Sensing what matters



Clinically tested*

Sleep monitor for detecting repetitive muscle jerking

The Emfit MM™, Movement Monitor is designed to detect faster, abnormal movements that may occur during the night, and will sound a high-frequency alert if these movements continue past a preset amount of time.

The Emfit MM, Movement Monitor consists of a film-like bed sensor and a control unit. The dynamic sensor pad is placed under the mattress and detects even the slightest movements of the person lying above it, including micro-movements caused by the heart beating. The sensor is attached to the control unit with an inter-connection cable. When faster, sustained movements continue longer than pre-set thresholds, a notification is triggered. System can also notice light movements, thus making it equally suitable for small children.

Emfit MM™, movement monitor is not a medical device and is not intended to be used for diagnosis or monitoring of disease or other conditions, or investigation of a physiological process.



Clinically proven, easy to install... and no false alarms

The top selling device of its kind and currently in use in over 60 countries worldwide, the Emfit MM is clinically proven technology that's easy to install. Sophisticated embedded software and patented sensor technology effectively prevent false alarms that can lead to user fatigue.

The Emfit MM has an audible alarm with adjustable volume for short-distance notifications. It can also interface with most nurse call systems, wireless transmitters and personal emergency phones with the included dry-contact output (NO and NC).

- **Clinically proven**
- **Patented, innovative sensor technology**
- **Detects micro-movements like heart beating from beneath the mattress**
- **Unobtrusive and non-body contact**
- **Durable, safe, and reliable**
- **No age or weight restrictions**
- **Easy to install**
- **Emits high-frequency alarm sound designed to carry to the caregiver**
- **Fast alert allows for a quick response time**
- **No false alarms**
- **Always on, always protected**
- **More sleep-filled nights and peace of mind**

Ultra-thin **bed sensor** is placed under the mattress



The Emfit Movement Monitor is designed as a nighttime motion detector for any type of sustained nighttime movement.

The control unit can be clipped to the bed using the included clip, mounted on the wall using the included fastening bracket, or placed at bedside. It is operated with 2 pcs AA size 1.5 V batteries. An optional, medical grade AC adapter is also available.

Clinical Trial Publications

Non-EEG seizure-detection systems and potential SUDEP prevention: State of the art

Anouk Van de Vel, Kris Cuppens, Bert Bonroy, Milica Milosevic, Katrien Jansen, Sabine Van Huffel, Bart Vanrumste, Lieven Lagae, Berten Ceulemans
Seizure 22 (2013) 345–355

Prospective Study of the Emfit Movement Monitor

Kate Van Poppel, Stephen P. Fulton, Amy McGregor, Michelle Ellis, Andrea Patters, and James Wheless
Journal of Child Neurology 28(11) 1434-1436, 2013.

Review Article – Sudden, Unexpected Death in Epilepsy

Orrin Devinsky
N Engl J Med 2011; 365:1801-11

Critical evaluation of four different seizure detection systems tested on one patient with focal and generalized tonic and clonic seizures.

Van de Vel A, Verhaert K, Ceulemans B.
Epilepsy Behav. 2014 Aug;37:91-4. doi: 10.1016/j.yebeh.2014.06.014. Epub 2014 Jul 7.
PMID: 25010322

Assessment of a quasi-piezoelectric mattress monitor as a detection system for generalized convulsions

Aditi P. Narechania, Irena I. Garic, Indranil Sen-Gupta, Mícheál P. Macken, Elizabeth E. Gerard, Stephan U. Schuele
Epilepsy & Behavior 28 (2013) 172–176

Testimonials



Emfit saved his life

Anyone considering EMFIT should definitely get one. We have had ours for a couple of months but never had a need until last night. Our son had his first tonic clonic in months at 2 am. His EMFIT alarmed and we were able to be with him, and put him on his side and help as his face was buried in his pillow. He had another with the paramedics present, but if it had not been for Emfit he may have stayed in status and we would not have known. Emfit saved his life.

Robert D.

The Emfit Monitor Has Changed Our Lives

The Emfit Monitor has changed our lives, I can't thank you enough! Our son is now able to be back sleeping in his own bedroom (and out of ours) now that we have a monitor to let us know if he has a seizure. I wish I could explain the relief we feel as parents to have something like this.

Jeanne Garske

We Can Sleep Easily for The First Time in 12 Years

This morning the monitor saved my son's life. We were alerted in the early hours of the morning that he was in the throes of a major status seizure. It was so early we would have never known he was having an issue - he could have easily aspirated as he had vomited during the seizure. We are so thankful that we were able to get one of your monitors and forever grateful. We can sleep easily for the first time in 12 years!

Anita Herold & Family – Beaverton, Oregon

The Device Is a God-send & Customer Service is Amazing

I cannot tell you what a relief this has been to me as a mother...but also to my son. He is a child and it worries him that we won't be there if he were to have a seizure. However, for all of our sanity we had to have a solution to sleepless nights. I have no qualms about recommending the Emfit Monitor. The device is a God-send and customer service is amazing.

Marie Smith – Laurel, Mississippi

We Credit the Emfit Monitor with Potentially Saving Our Daughter's Life

We purchased an Emfit Movement Monitor on our own from one of your overseas suppliers a little over a year ago. Our daughter has seizures that have never stopped on their own and almost always at night (stopped by Diastat). They are, thankfully, fairly rare, occurring perhaps once every couple of months. On at least two occasions, when we didn't catch the seizure with the baby monitor (for various reasons), we credit the Emfit monitor with potentially saving our daughter's life.

Jeremy Hilton – Montgomery, Alabama

His life has most definitely been saved a few times

I am so grateful to you all at Emfit. I have three of your alarms and all of my boys can sleep soundly now. My eldest says he knows he can sleep safely and confident in knowing he's safe! This in itself has brought a warm tingle to me knowing that he feels like this...his life they said, has most definitely been saved a few times. Thank you all for such a life-saving alarm.

Lee-Anne M.

It gives us peace of mind

We recently purchased a monitor for our 15-year-old daughter who has tonic clonic generalized seizures. Last night at 8pm the monitor went off while she was sitting on her bed working on homework. I would have never known she was having a seizure if the alarm had not gone off. Then, the alarm went off again at 5 am and we were alerted to another seizure. As a 15-year-old, she does not want to have to keep her door open all the time or have mom or dad sleep in her room. So thankful for our monitor. It gives us peace of mind that we will be alerted if she is having a seizure, while allowing her to have her privacy and not have to feel like we are hovering over her all the time. Thank you Emfit.

Laura R.

www.emfit.com

EMFIT

USA

Epi USA

415 N. Guadalupe St., #403
San Marcos, TX 78666
Phone: (512) 667-7507
christine.ocean@emfit.com
www.epiusa.net

EUROPEAN UNION

Emfit Ltd

Konttisentie 8 B
40800 Vaajakoski
FINLAND
Phone: +358 20 778 0870
info@emfit.com

Designed and manufactured in Finland